




# LSM1301 General Biology

## *Laboratory 2: DNA Extraction, PCR, and Gel Electrophoresis*

*Semester 1, AY2025/2026*



# STARTING TASKS



- Use alcohol and Kimwipe to clean your goggles
- Goggles and gloves needed for today's lab work

# LEARNING OBJECTIVES

1. Isolate DNA from banana and spinach
2. Explain how the cell is lysed and DNA is released from cells
3. Conduct PCR and explain the principle of PCR
4. Conduct agarose gel electrophoresis and explain its principle
5. Interpret experimental data and to report data effectively

# ORDER OF TASKS

1. Gel Electrophoresis of PCR samples (*45 mins to 1 hour*)
2. DNA Extraction from banana and spinach (*20 mins*)
3. View and analyse gel electrophoresis results (*5 to 10 mins*)

# Previously...

Each pair of students prepared a PCR sample containing the following:

- 10  $\mu$ l PCR master mix
- 2  $\mu$ l Forward primer (*note down set number*)
- 2  $\mu$ l Reverse primer (*note down set number*)
- 2  $\mu$ l DNA template (plasmid)
- 4  $\mu$ l Distilled water

*Subjected to the following PCR conditions:*


STEP	TEMP.	TIME
Initial Denaturation	95°C	3 min
Denaturation, Annealing, and Extension (35 cycles)	95°C	30 secs
	55°C	30 secs
	72°C	30 secs
Final Extension	72°C	5 min
Hold	15°C	5 min

# TODAY: Gel Electrophoresis

## Purpose: To find out if the PCR was successful

- Gel electrophoresis separates molecules by size and charge
- Agarose gel contains staining agent called SYBR Safe
- DNA molecules visualised under UV light after electrophoresis
- DNA standards (“ladder” or “marker”) as reference sizes

GeneRuler 100 bp Plus DNA Ladder




0.5 µg/lane, 8 cm length gel,  
1X TBE, 5 V/cm, 1 h

5% polyacrylamide

0.5 µg/lane, 20 cm length gel,  
1X TAE, 8 V/cm, 3 h

# TASK 1: Gel Electrophoresis

- TA will demonstrate sample loading using 15  $\mu$ l DNA ladder
- Students will load 15  $\mu$ l of their own sample
- Run electrophoresis at 120V for 40 minutes
- Visualise gel under UV light
- Take a photo; use it to answer Q2 of assignment



## TASK 2: DNA Extraction



Starting materials:  
quarter slice of banana; one spinach leaf



Place in mortar



Add a few drops of salt water; grind into paste  
with pestle



Add remaining salt water; use pestle to gently  
mix water + paste

*Individual Work: each student should prepare one banana  
extraction and one spinach extraction.*



## TASK 2: DNA Extraction

Transfer contents to tube containing liquid soap



Cover; invert gently to mix for 2 to 3 mins



Tilt tube; gently add **ice-cold** alcohol down one side of the tube




Cover; let tube rest for 5 mins (don't shake!)



DNA will float to the top of solution



To store DNA: use glass rod to twirl DNA, place in 1 mL **room temp.** alcohol



# TASK 3: How to interpret/analyse gel photo

## 1. Label lanes for easy reference

*Legend:*

*1: DNA marker*

*2: Sample 1*

*3: Sample 2*

*4: ...*


*5: ...*



# How to interpret/analyse gel photo

## 2. Identify and label DNA marker sizes

GeneRuler 100 bp Plus DNA Ladder



1000 bp →

500 bp →

**Important: Label the relevant marker sizes in the gel photo, not just the two examples shown in this photo.**




# How to interpret/analyse gel photo

1 2 3 4 5 6

3. Estimate size of sample bands in all lanes

*e.g. The bright band in lane 3 is between 200 to 300 bp*

**Note:** It is not necessary to draw the dotted lines around the bands, or lines aligning to the marker sizes. These lines are for demonstration only to show you how to estimate the band sizes. Avoid devising in-between sizes, such as “250 bp”, which is not one of the reference sizes in the DNA ladder.



# CITATIONS: How To?

FYI



For certain lab assignment questions, you can cite information you obtain from scientific articles:

- You can use Google scholar (<https://scholar.google.com/>) or the NUS library search engines (<https://nus.edu.sg/nuslibraries/databases-search>) to do this
- Remember to rephrase the article's text into your own words and to cite the paper properly (below your answer).
- You may use APA citation format :
  - Author Last name, First initial. Middle initial. (Year Published). Title of article. *Title of Periodical*, Volume(Issue), page range.
- Example
  - Ross, R. P., Morgan, S., & Hill, C. (2002). Preservation and fermentation: past, present and future. *International journal of food microbiology*, 79(1-2), 3-16.